


# No. 220 High Speed Air Turbine Laboratory Handpiece System



*TOTALLY NEW Vibration Free and Super Quiet!*

**Powerful, Vibration Free, Super Quiet, and a Smooth Cutting Performance in a Dependable and Economical Package!**



- **Decades of Success**

Time-proven workhorse air turbine in use in thousands of laboratories around the world

- **High-Powered Cutting Action**

350,000+ RPM provides powerful cutting action with standard FG burs

- **Ergonomic Design**

Small, light-weight aluminum design is easy to handle with less hand fatigue

- **Rugged Turbine Cartridge**

Durable friction Cartridge/Chuck system features long term reliable performance with gross trimming and coarse cutting tasks

- **Simplified Maintenance Requirements**

New Insertion and Removal Tool Sets simplify Cartridge and Metal Chuck maintenance requirements

- **Small Size**

4.25" Long x .68" Diameter x 2.2 oz.

- **Ordering Information**



No. 220 Complete System  
Also Includes Filter/Regulator,  
Lubricator, Tubing and Oil

| <b><i>Item #</i></b> | <b><i>Description</i></b>   |
|----------------------|---|
| 01000                | No. 220 Air Turbine Complete System (Includes Handpiece, Regulator, Lubricator, Foot Control, Tubing and Oil) |
| 01108                | No. 220 Air Turbine Handpiece Only  |
| 01116                | Replacement Turbine (Cartridge) for No. 220   |
| 01120                | Replacement Metal Chuck for No. 220 (and MaxAir)  |
| 01122                | Installation (Compression) Tool for Metal Chuck   |
| 01121                | Standard Removal Tool for Metal Chuck   |
| 01125                | Pro Removal Tool Set for Metal Chuck (Hardened Screw + Steel Punch Set)                                       |
| 01109                | Turbine Oil, Low Viscosity, 1 Oz. Bottle  |
| 01127                | Turbine Oil, Low Viscosity, 1 Qt. Bottle  |

## **Buffalo Dental Manufacturing Co., Inc.,**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

**Tel:** 516.496.7200 or 800.828.0203 **Fax:** 516.496.7751

**[www.buffalodental.com](http://www.buffalodental.com)**

S01000 1109

***Please order today from  
your favorite Dealer!***