

The Dynamic Duo



***Two of our most popular units
combined into one complete system!***

*Patents Pending

***A high-powered, reliable Electric
Laboratory Handpiece System
and an effective, time-saving
knife all at your fingertips!***



Features:

- Four Multi-Speed options (14k, 22k, 27k, 35k rpm) allows for ideal speed selection for most cutting needs.
- High-Torque motor provides smooth operation and powerful performance with less operator fatigue.
- Fully adjustable temperature control for precise cutting action. Cutting tip heats ultra-fast to save time and eliminate waiting.
- A & B selection switch allows for the use of two different Thermal Pens to reduce time and increase trimming productivity.
- Lightweight pen features an ergonomic design and sharp, permanent tip for improved comfort and control.

Ordering Information:

| Item# | Description | |
|--------------|---|--|
| 38760: | Dynamic Duo Complete System | 80510: Replacement Pen w/ Concave Tip |
| 37590-14: | Replacement Metal Chuck for Gray Handpiece | 80520: Replacement Pen w/ Knife Edge Tip |
| 37590-37: | Replacement Carbon Brushes for Gray Handpiece (2) | 80530: Replacement Pen w/ Waxing Tip |

Distributed by:

Buffalo Dental Manufacturing Co., Inc.
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

*Please order today from
your favorite dealer!*