

MOLLOPLAST-B®

Permanent Soft Relining Material




MOLLOPLAST-B® is an industry-leading heat-processed silicone material that stays soft permanently and features decades of documented success in denture reline cases. The material bonds firmly to new or existing acrylic dentures, and stays elastic and bacteria free for years



Features/Benefits:

- Single component silicone
- Heat-curing
- Easy to use with dental flasks: can be polymerized either in water bath or microwave
- Strong bonding to any denture base, to new and old dentures alike
- Time saving feature: can be polymerized simultaneously with acrylates
- Resistant to the influences of the oral environment
- Contains no substances of swelling and no plasticizers
- Non-irritant, tissue-compatible, clinically tested and scientifically documented
- MOLLOPLAST-B® relinings will remain soft and elastic even after years of wearing



Ordering Information:

<i>Item #</i>	<i>Description</i>
62300	MOLLOPLAST-B®, 45 gm Reg. Size
62310	MOLLOPLAST-B®, 170 gm Lab Size
73000	Primo, 15 ml Bottle and Brush

MOLLOPLAST-B®, the industry leader in the field of soft liner materials.

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751

www.buffalodental.com

Please order today
from your favorite dealer!