

# Buffalo Dental Waxing Unit



## ***Compact and Economical 3-Compartment Wax Heating Unit with Adjustable Temperature Control***

### ***Features/Benefits***

- On/Off switch located on power cord to prevent wax seepage into switch mechanism
- Controlled temperatures from 140° to 210° F for convenient operation (adjustment screw on back left)
- Provides three 11/16" deep wax compartments (one is 2-1/2" x 2-7/8", and two are 1-5/8" x 1-3/8")
- Attractive design with finished metal sides, handles and helpful heat indicator light on front panel



### ***Ordering Information***

<u>Item No.</u>	<u>Description</u>
86000	Waxing Unit, 120V AC
86010	Waxing Unit, 220V AC

**Buffalo Dental Manufacturing Co., Inc.,**  
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA  
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751  
[www.buffalodental.com](http://www.buffalodental.com)

S86000 0706

*Please order today from  
your favorite Dealer!*