

# ThermaPlier™

## Hot Tipped Thermal Pliers



*For accentuating and making retentive undercuts on orthodontic appliances and mouthguards*

Reheats the critical retention area allowing emphasis without causing adjacent plastic to distort.

- Adjust the temperature to a setting which allows the creation of dimples in the material you are working with but which does not burn the material.
- Use the plier to make retentive undercuts on the orthodontic appliances and mouthguards as per figure 1 & 2.

\*Patent Pending.

*Made in USA!*



Figure 1



Figure 2

### • Ordering Information

<u>Item #</u>	<u>Description</u>
80610	ThermaPlier 110 V

Distributed by:

Buffalo Dental Manufacturing Co., Inc.  
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA  
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751  
[www.buffalodental.com](http://www.buffalodental.com)

S18100 0209

*Please order today from  
your favorite Dealer!*