

Econo-Vac™ Vacuum Forming System



Classic High Performance Vacuum Forming System in Steady and Successful Service in Thousands of Dental Offices and Labs

- **High Performance System**

Dual-stage electric vacuum motor provides powerful and immediate suction... and highly accurate adaptations- no compressed air connections required!

- **Accurate Frame and Alignment System**

Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

- **Heavy Duty Construction**

Solid construction (7 1/2" W x 10 1/2" D x 15" H, 32 lbs.) for maximum durability, long life and virtually no maintenance!

- **Convenient and Safe Design**


Heater hood rotates behind the center post for safety and convenience, and the compound frame closure latch helps keep the components cool to the touch during repeated use

- **Wide Range of Available Accessories**

Please see the complete line of Buffalo Dental vacuum forming equipment and supplies including a wide range of Vacuum Forming Sheets, Block-out Compound, Silicone Spray, Tray Handles, Knives and Burs, the ThermaKnife™ Thermal Cutting Knife, and PolyShoK™ Sports Mouthguard Materials



- **Ordering Information**

Item #	Description
80186	Econo-Vac Vacuum Forming System Complete, 120V AC, incl. a reference assortment of vacuum forming materials
80187	Econo-Vac Vacuum Forming System Complete, 220V AC,  incl. a reference assortment of vacuum forming materials

For laboratory use only.

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 **Fax:** 516.496.7751
www.buffalodental.com

580186 0708

**Please order today from
your favorite Dealer!**