

Stalite™ Special Solution



Unique solution helps make harder, smoother, and more accurate stone models & dies

- Use instead of water with gypsum to create models and dies that are consistently harder, stronger, smoother and more abrasion resistant
- Creates harder models with more flexural strength to resist breaking, chipping and abrasion
- Saves time by promoting consistent, thin mixes with less bubble formation
- Smoother final surfaces help improve model and restoration accuracy
- Improved freeze-stable formula
- Can be used with all gypsums
- No significant setting expansion

Item #	Description
80151	Stalite Special Solution, 1 Gallon Bottle



Distributed by:

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

580151 0706

**Please order today from
your favorite Dealer!**