

BUFF-TUFF **Mouthguard Laminate**



Introducing our ***NEW BUFF-TUFF***
Mouthguard Laminate!

NO BPA!

NO DIOXIN!

NO PVC OR INTERNAL PVC BRACES!

CROSS-GRAINED FOR STRENGTH

AND CONSTANT THICKNESS OF

FINISHED MOUTHGUARD!

DOES NOT PRODUCE OFFENSIVE SMOKE OR TOXIC FUMES!

DOES NOT SEPARATE!



BUFF-TUFF Mouthguard Laminate

*Patents Pending.

BUFF-TUFF Mouthguard Laminate features Cross-Graining for extra strength and constant thickness of finished mouthguard. Injection molded mouthguards cannot be cross grained.

-Designed Specifically For Athletic Use-

Ordering Information:

| <i>Item#</i> | <i>Description</i> |
|--------------|--|
| 70000 | <i>BUFF-TUFF</i> .160" Clear Mouthguard Material, 10 Square Sheets |
| 70010 | <i>BUFF-TUFF</i> .160" Clear Mouthguard Material, 100 Square Sheets |
| 70100 | <i>BUFF-TUFF</i> .160" Clear Mouthguard Material, 10 Discs |
| 70110 | <i>BUFF-TUFF</i> .160" Clear Mouthguard Material, 100 Discs |

****Also available in Blue, Black, Orange, Yellow, Red, and White. Call for more information.***

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203
Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com