

Magic Wand™

Electric Lab Handpiece System



Versatile and Compact 35,000 RPM System Ideal for Medium Duty and Intermittent Lab Trimming Needs

- Compact electric handpiece system features fully variable cutting speeds from 0 to 35,000 RPM
- Premium black electric handpiece provides smooth, quiet and comfortable operation with a history of long term durability
- Reliable, no-frills, made-in-USA Console
- Fully variable speed conveniently controlled by a heavy-duty, all-metal Variable Speed AC Foot Control
- The convenience and benefits of an electric handpiece without the added extras of higher priced systems



• Sizes

Handpiece: 6 in. (153 mm) Long x 1 1/8 in. (29 mm) Major Diameter, 5/8 in. (16 mm) Nose Diameter x 9.5 oz. (270 gm) Weight

Console: 4.25 in. (108 mm) Width x 5.25 in. (133 mm) Depth x 3 1/8 in. (80 mm) Height, 2 lb. 7 oz. (1.1 kg) Weight

• Ordering Information

<i>Item #</i>	<i>Description</i>
38500	Magic Wand Electric Handpiece System Complete, 120V AC, including Silver Control Console, Handpiece w/HP Collet, Foot Control, Spare Square Carbon Brushes- Set of 2, Handpiece Holder and Collet Wrenches
38600	Magic Wand Electric Handpiece System Complete, 220V AC, including Silver Control Console, Handpiece w/HP Collet, Foot Control, Spare Square Carbon Brushes- Set of 2, Handpiece Holder and Collet Wrenches
38025	Black Signature Electric Lab Handpiece Only
38025CN	Replacement HP Collet Chuck for Handpiece, 3/32" Size
38025BS	Replacement Square Carbon Brushes for Black Signature Electric Lab Handpiece, Set of 2

Note: Use caution and wear personal protective eyewear and equipment when operating the Magic Wand system

Buffalo Dental Manufacturing Co., Inc.,
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | **Tel:** 516-496-7200 or 800-828-0203
Fax: 516-496-7751 | **Email:** Sales@mailbdm.com | www.buffalodental.com
S38500 0418