

# Model Arch Trimmer



**Premium High-Torque Trimmer Provides Fast and Quiet Cutting for the Most Demanding Model Arch Dry Trimming Applications!**



- **Fast and Easy Dry Trimming**

High-torque, 5400 rpm, 1/8 HP heavy-duty motor promotes quiet, fast, and nearly-effortless cutting

- **Improved Dust Control**

Enlarged 7" x 7" work table with built-in dust control ducting and 1" OD exhaust port on back maintain quiet and effective dust control- ready for attachment to your existing dust collector

- **Rapid Tool Changes**

Quick-release work table and 1/4" collet allow convenient and fast changing of 1" maximum OD cutting tools

- **Improved Safety**

Strong 4" x 6" clear polycarbonate safety shield helps protect the operator from debris, and the interlocking on-off switch automatically resets to off in case of power interruption

- **Buffalo Tough Cutting Instruments**

System ships complete with a fast cutting Buffalo **Tough** 1/2" diamond-cut tungsten carbide bur and a 3/4" coarse-grit arbor band set

- **Compact Size**

8 in. (20.4 cm) Width x 10 in. (25.4 cm) Depth x 8-1/8 in. (20.6 cm) Height x 14 lb. (6.4 kg) Weight

- **Complete Range of Products**

Please see our complete range of model trimming equipment and supplies including Wet and Dry Model Trimmers, Trimming Wheels and Belts, Knives and Burs, Safety Shields, and Disposable Plaster Traps.

- **Ordering Information**



**Made in USA!**

<b>Item #</b>	<b>Description</b>
64000	Model Arch Trimmer Complete System, 120V AC, 1.2 A, including Safety Shield, Wrench Set, 1/2" Diameter Diamond-Cut Tungsten Carbide Bur (#19400), 3/4" Arbor Band Mandrel (#02925), and Coarse Grit Zirconia Arbor Bands, Pkg. of 10 (#02700)

*Notes: For laboratory use only. Use extreme caution and wear personal protective eyewear, gloves and equipment when operating a Model Arch Trimmer. Avoid personal contact with the spinning bur, and always unplug the Trimmer when changing burs or cleaning the housing.*

**Buffalo Dental Manufacturing Co., Inc.,**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | **Tel:** 516-496-7200 or 800-828-0203

**Fax:** 516-496-7751 | **Email:** [Sales@mailbdm.com](mailto:Sales@mailbdm.com) | [www.buffalodental.com](http://www.buffalodental.com)

S64000 0418