

# Vibrators



*Our Premium Vibrators come in a variety of options. Ranging from small, medium and large, Buffalo Dental Vibrators are powerful and effective!*

## No. 200 Extra - Heavy Duty Vibrator

Our Largest, Most Powerful Vibrator for Multiple Flasks and Dense Castings

Large 10 ½" long x 7 ½" wide rubber-covered work platform for convenience and spill control. Improved electronics and fully variable speed control for easy and powerful operation. Solid, beige enameled housing provides stability and eases clean-up tasks. 4 ½" high. 19 lbs.

**84500** No. 200 Vibrator, 120 V AC



## No. 2 Heavy Duty Vibrator

Medium Size with High Powered Vibration

Ribbed 6" diameter rubber-covered platform designed to hold larger flasks. Features a three-speed control switch for variable operation. Beige enameled housing (11" wide x 4 ½" deep x 6 ½" high, 13 lbs.) for stability and performance.

**84400** No. 2 Vibrator, 120 V AC

## No. 1A Vibrator

Our Most Popular Vibrator- At a Surprisingly Affordable Price!

Improved 3-speed electronics and strong, dependable, performance make this our most popular small vibrator. Durable 4" diameter rubber cover and suction cup feet improve stability and ensure quiet operation. Small size (4 ½" long x 4 ½" wide x 4 ½" high, 6 lbs.) fits conveniently in office labs.

**84350** No. 1A Vibrator, 120 V AC

**84324** SplashCap™ Splash Guard for 1A and 1B Vibrators



## No. 1B Vibrator

Small Vibrator with Variable Speed Foot Control

Identical size to our 1A Vibrator, but features foot pedal operation instead of a rotary switch. Pedal offers convenient hands-free operation and fully variable speed control.

**84360** No. 1B Vibrator, 120 V AC

**Buffalo Dental Manufacturing Co., Inc.,**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516-496-7200 or 800-828-0203

Fax: 516-496-7751 | Email: [Sales@mailbdm.com](mailto:Sales@mailbdm.com) | [www.buffalodental.com](http://www.buffalodental.com)