


# VACUUM FORMING EQUIPMENT & MATERIALS



Buffalo Dental Manufacturing Company, Inc.  
[www.buffalodental.com](http://www.buffalodental.com)



## The Dynamic Duo™

A time-saving knife with a hot-tipped blade that promotes fast and easy trimming of ALL vacuum formed materials. Adjustable temperature control for precise heat adjustment providing the operator with immediate heat for fast cutting and trimming. Ultra-fast heating circuit gets tip to ideal temperature in seconds, saving time. Lightweight pen with sharp, permanent tip promotes improved comfort and control. We have even added an EXTRA OUTLET for an additional Pen. Unit Comes complete with 1 Knife edged Tip Pen.

### Item # Description

- 38760 Buffalo Dynamic Duo V-35/Therma-Knife Combo Unit Complete, 120 V AC  
 38765 Buffalo Dynamic Duo V-35/Therma-Knife Combo Unit Complete, 220 V AC



The Dynamic Duo™

## V-35 MicroMotor – Electric Handpiece System

The V35 electric handpiece system features a no-frills, multi-speed control console teamed with our new Gray Handpiece. It represents a heavy-duty, maximum 35,000 RPM, lab-only handpiece system for use with HP-sized burs that offers high durability, high torque, and economical pricing.

### Item # Description

- 37570 V35 Gray Electric Handpiece System, 120V AC (V35 Console, Gray Electric Handpiece, Foot Control and Handpiece Rest)  
 37575 V35 Gray Electric Handpiece System, 220V AC (V35 Console, Gray Electric Handpiece, Foot Control and Handpiece Rest)



V-35

## ThermaKnife™ Hot Tipped Tray Trimming Knife

Effective and Time-Saving Knife For Fast and Easy Trimming of Bleaching Trays, Night Guards, Mouthguards, and Other Thermoplastic Materials.

- Fully adjustable temperature control for precise cutting action
- Cutting tip heats ultra-fast to save time and eliminate waiting
- Lightweight pen features an ergonomic design and sharp, permanent tip for improved comfort and control Knife-Edge Tip
- A & B selection switch allows for the use of two different Thermal Pens to reduce time and increase trimming productivity
- Fast and easy pen replacement (via snap connector at back of pen) for added convenience Concave Tip
- Convenient pen rest provides added safety when unit not in use
- Small portable size: Console 5.75" W x 5" D x 4"H (incl. Rest) Pen: 6" L x .75" W (at foam grip,) 2.5 lb. total console/pen weight



ThermaKnife w/ Knife Edge Tip

### Item # Description

- 80500 ThermaKnife Intro Kit, 120V AC  
 (Includes Knife Edge Pen)
- 80510 Replacement Pen with Concave Tip  
 80520 Replacement Pen with Knife Edge Tip  
 80530 Replacement Pen with Waxing Tip

## ThermaPlier™ Hot Tipped Thermal Pliers

For accentuating and making retentive undercuts on orthodontic appliances and mouthguards. Reheats the critical retention area allowing emphasis without causing adjacent plastic to distort.

- Adjust the temperature to a setting which allows the creation of dimples in the material you are working with but which does not burn the material.
- Use the plier to make retentive undercuts on the orthodontic appliances and mouthguards as per figure 1 & 2.



Figure 1



Figure 2



ThermaPlier™

### Item # Description

- 80610 ThermaPlier 110V

### Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA  
 Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax/orders: (516) 213-6610  
 www.buffalodental.com

### Accu-Vac™ Precision Vacuum Forming System

An extremely powerful, accurate, and quiet vacuum forming system utilizing compressed air-powered suction.

- Utilizes ultra-high powered venturi suction (requires 60 to 100 psi external compressed air supply) to create highly accurate vacuum adaptations
- Gasket-free frame accepts 5" diameter round and 5" x 5" square materials for improved versatility

- Unique table uses either standard flat plate or included steel bead basket for multiple model positioning options
- Single On/Off heat control switch and automatic vacuum actuation switch promote easy operation
- Convenient small size (7 1/2" W x 10 1/2" D x 10" H x 16 lbs.) enhances portability
- System includes Accu-Vac unit, Flat Plate Table, Steel Beads, Bead Basket, and sample range of Vacuum Forming Materials

Item #	Description
80175	Accu-Vac System, 120 V AC
80176	Accu-Vac System, 220 V AC



Accu-Vac

### Sta-Vac II™ Versatile Vacuum Forming System

Our improved classic high performance electric vacuum former that now allows both traditional flat table and highly accurate bead basket adaptations.

- New ConverTable™ combination vacuum table/bead basket can be used in standard flat table mode or reversed and used in bead basket mode for improved versatility, vacuum flow, and accuracy
- Dual-stage electric vacuum motor provides powerful suction and highly accurate adaptations- no compressed air connections required!
- Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table
- Solid construction (7" W x 16" D x 11" H x 32 lbs.) for maximum durability and long life
- Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation

- All systems include a sample range of vacuum forming materials

Item #	Description
80193	Sta-Vac II Vacuum Forming System, 120 V AC, w/ Two ConverTables, Steel Beads and Plastic Work Tray
80194	Sta-Vac II Vacuum Forming System, 220 V AC, w/ Two ConverTables, Steel Beads and Plastic Work Tray

Also available is the Sta-Vac, with all internal components identical to the Sta-Vac II, but with only one ConverTable and no steel beads nor work tray, designed for use in original flat model table mode only.

Item #	Description
80191	Sta-Vac System, 220 V AC
80190	Sta-Vac System, 120 V AC



Sta-Vac


### Econo-Vac™ Vacuum Former

Premium vacuum former features power and performance in a standard, upright design.

- Dual-stage vacuum motor provides powerful suction and accurate adaptations
- Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials

- Square center post assures smooth operation and guarantees perfect alignment of table and frame
- Solid construction (7" W x 10 1/2" D x 15" H x 32 lbs.) for high durability and long life
- All systems include a sample range of vacuum forming materials

Item #	Description
80186	Econo-Vac System, 120 V AC
80187	Econo-Vac System, 220 V AC



Econo-Vac

### Tray-Vac™ Vacuum Former

A high performance, no-frills, economical, and reliable small vacuum former for the dental office.

- High-velocity air flow system promotes rapid, highly accurate adaptations with 5" x 5" sheet materials
- Unique quick-lock frame and gasket-free design promote convenient, reliable operation

- Lightweight, small machine (7 1/2" W x 9 1/2" D x 15" H x 17 lbs.) improves portability
- All systems include a sample range of vacuum forming materials

Item #	Description
80165	Tray-Vac, 120 V AC
80168	Tray-Vac, 220 V AC



Tray-Vac



## Sta-Lube™ Spray

A convenient silicone lubricant spray that prevents materials from sticking to models, vacuum table and frame.

- Environmentally safe material assures easy clean-up of flasks and trays
- Ideal separating medium for dies when making copings

**Item # Description**  
80015 Sta-Lube Silicone Spray, 16 oz. Can



Sta-Lube

## Alcohol Torch

An economical torch that produces a needle-point flame with only a gentle squeeze of the plastic bottle.

- 5" high torch uses denatured alcohol fuel for cost savings and convenience
- Base of bottle weighted with marbles to prevent tipping and improve safety
- Ideal for applying a rolled edge to vacuum formed adaptations (flame edge and rub with thumb)
- Can also be used as an alcohol lamp

**Item # Description**  
82440 Alcohol Torch  
57560-1D Wicks for Alcohol Torch, Pkg. of 12



Alcohol Torch

## No. 55 Plate Shears

Heavy duty scissor-style instruments for multiple cutting and trimming applications.

- No. 61 plate shears excellent for trimming vacuum forming materials

**Item # Description**  
78250 No. 55 Plate Shears



#55

## Arbor Band Mandrels

High quality HP (3/32" diameter) mandrels with adjustable -diameter rubber heads are designed to firmly hold arbor bands in place during use.

- Replacement rubber heads available for renewing older mandrels

**Item # Description**  
02771 3/8" Medium Arbor Bands , pkg 100  
02920 3/8" Complete Mandrel, HP Sized Shank



## Burs, Miscellaneous

High quality burs for miscellaneous applications.

- Cutting and Finishing Bur utilizes hardened tool steel for fast and accurate cutting to 15,000 rpm. Bur is ideal for trimming custom trays, base plates, and tough thermoplastic sheet materials.

- Plastic & Stone Trimming Burs designed for fast and smooth trimming of stone and epoxy dies, models, acrylics, and sheet resin materials. Burs constructed of carbide steel for long life. Popular taper shape with non-clogging flutes.

**Item # Description**  
20000 Cutting & Finishing Bur for Thermoplastic Materials, HP



Cutting & Finishing Bur

## PRKIT

**Item # Description**  
18100 Cone 82T carbide bur  
18010 Taper 84T carbide bur  
02781 3/8" Coarse arbor bands (100)  
02920 3/8" Complete arbor band mandrel

## Block-Out Compound

Used for blocking out remaining teeth and undercuts on models when making appliances.

- Reusable, heat resistant material resists running and sticking

- Can also be rolled around the base of models to promote easy periphery trimming

**Item # Description**  
05240 Block-Out Material, 1 lb.

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax/orders: (516) 213-6610  
www.buffalodental.com



## Vacuum Forming Materials

High quality 5" x 5" sheet thermoplastic materials for creating a wide range of vacuum formed appliances and trays.

- Pure grade resin sheets are latex-free and uniform in size, cut and thickness for improved performance.
- Resist bubbling, crazing and fracturing to promote superior adaptations every time.
- Clear soft Ethyl Vinyl Acetate (EVA) sheets designed to bond together for build-up and lamination of custom mouthguards.

### \*NEW\* Clear Tray Material

Buffalo Clear Material offers the best model adaptation and excellent rigidity. The result will be a clear, highly polished surface with no color dyes or bluing agents. The Extra Wide size stays secure in all frames.

#### .020" (0.5mm) Clear Tray Material

Item #	Description			
81660	Pkg. of 50	81665	Pkg. of 200	81670 Pkg. of 625

#### .040" (1.0mm) Clear Tray Material

Item #	Description		
30990	Pkg. of 25	30995	Pkg. of 625

#### .060" (1.5mm) Clear Tray Material

Item #	Description		
31000	Pkg. of 25	31010	Pkg. of 625

#### .080" (2.0mm) Clear Tray Material

Item #	Description		
31090	Pkg. of 25	31100	Pkg. of 625

#### .100" (2.5mm) Clear Tray Material

Item #	Description		
31110	Pkg. of 25	31120	Pkg. of 625

#### .020" Coping and Temporary C&B Material

Cloudy, stiff Polypropylene sheets excellent for use when fabricating temporary crowns, temporary bridges, and copings. Sized for improved vacuum is easier to use and promotes improved vacuum seal. Will not adhere to acrylics.

Item #	Description			
31720	Pkg. of 50	31725	Pkg. of 200	31730 Pkg. of 625

#### .060" Pink Base Plate Material

Stiff pink Polystyrene sheets excellent for fabrication of thin base plates.

Item #	Description			
04610	Pkg. of 25	04615	Pkg. of 200	04620 Pkg. of 625

#### .080" Pink Tray Material

Stiff pink Polystyrene sheets excellent for fabrication of medium and thick base plates and custom trays.

Item #	Description			
82610	Pkg. of 25	82615	Pkg. of 200	82620 Pkg. of 625

#### .100" Pink Tray Material

Stiff pink Polystyrene sheets excellent for fabrication of thick base plates and custom trays.

Item #	Description			
82650	Pkg. of 25	82655	Pkg. of 200	82660 Pkg. of 625

#### .125" Blue Tray Material

Stiff blue Polystyrene sheets ideal for fabrication of thinner custom trays.

Item #	Description			
82690	Pkg. of 25	82695	Pkg. of 200	82700 Pkg. of 625

#### .150" Blue Tray Material

Stiff blue Polystyrene sheets ideal for fabrication of thick custom trays.

Item #	Description			
82730	Pkg. of 25	82735	Pkg. of 200	82740 Pkg. of 625

#### .040" Clear Bleaching Tray Material

Clear, soft EVA Sheets excellent for use when fabricating bleaching trays. Sized for improved vacuum is easier to use and promotes improved vacuum seal.

Item #	Description					
62820	Pkg. of 12	62803	Pkg. of 25	62830	Pkg. of 100	62840 Pkg. of 300

#### .080" Clear Thin Mouthguard Material

Clear, soft EVA Sheets excellent for fabrication of thin mouthguards for non-contact sports, bruxism splints, soft splints, and custom fluoride trays. Material bonds to itself so it can be built up in layers for creation of Athletic Mouthguards. Contains *NO PVC*.

Item #	Description				
62804	Pkg. of 12	62806	Pkg. of 25	62802 Pkg. of 100	62800 Pkg. of 300

#### .118" Clear Mouthguard Material

Clear, soft EVA Sheets excellent for fabrication of light mouthguards, soft splints, and custom fluoride trays. Material bonds to itself so it can be built up in layers for creation of Athletic Mouthguards.

Item #	Description				
62815	Pkg. of 12	62816	Pkg. of 25	62805 Pkg. of 100	62810 Pkg. of 300

#### .160" Clear Mouthguard Material

Thick, clear, resilient EVA sheets for fabrication of thick mouthguards and resilient bruxism splints. Material bonds to thin EVA mouthguard materials so adaptations can be built up in layers for creation of Athletic Mouthguards.

Item #	Description				
62700	Pkg. of 12	62725	Pkg. of 25	62705 Pkg. of 100	62710 Pkg. of 300

#### .160" Sta-Guard™ Clear Mouthguard Material

Thick, clear, resilient EVA classic-brand sheets excellent for fabrication of thick mouthguards and resilient bruxism splints. Material bonds to thin EVA mouthguard materials so adaptations can be built up in layers for creation of Athletic Mouthguards.

Item #	Description			
62750	Pkg. of 12	62755	Pkg. of 100	62760 Pkg. of 300

### Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com



## Thick resilient 5" X 5" EVA Colored Mouthguard Sheets, Excellent Protection for All Sports!

Buffalo Dental Mouthguard & Tray material feature a "Total Smooth" surface for comfort and appearance! Contains *NO BPA, NO PVC, NO LATEX, or PVC INSERTS!* Also features NO inclusions or glitter which can be swallowed! This material is used for Athletic Mouthguards.

**.040" EVA Mouthguard Material** 62610 PINK, Pkg. of 12

Item #	Description
72560	RED, Pkg. of 12
72570	RED, Pkg. of 100
72580	RED, Pkg. of 300
72590	BLUE, Pkg. of 12
72592	BLUE, Pkg. of 100
72593	BLUE, Pkg. of 300
72595	ORANGE, Pkg. of 12
72596	ORANGE, Pkg. of 100
72597	ORANGE, Pkg. of 300
72585	BLACK, Pkg. of 12
72586	BLACK, Pkg. of 100
72587	BLACK, Pkg. of 300
72581	WHITE, Pkg. of 12
72582	WHITE, Pkg. of 100
72583	WHITE, Pkg. of 300
72476	YELLOW, Pkg. of 12
72477	YELLOW, Pkg. of 100
72478	YELLOW, Pkg. of 300

**.040" EVA Mouthguard Material Disc (120mm)**

Item #	Description
62827	CLEAR, Pkg. of 12

**.040" EVA Mouthguard Material Disc (125mm)**

Item #	Description
62827	CLEAR, Pkg. of 12
62822	CLEAR, Pkg. of 100
62825	CLEAR, Pkg. of 300
72670	RED, Pkg. of 12
72680	RED, Pkg. of 100
72690	RED, Pkg. of 300
72691	BLACK, Pkg. of 12
72692	BLACK, Pkg. of 100
72693	BLACK, Pkg. of 300



MADE IN USA

**.080" EVA Mouthguard Material**

Item #	Description
62788	BLACK, Pkg. of 12
62795	BLACK, Pkg. of 300
62789	WHITE, Pkg. of 12
62799	WHITE, Pkg. of 300
62784	ORANGE, Pkg. of 12
62798	ORANGE, Pkg. of 300
62792	YELLOW, Pkg. of 12
62794	YELLOW, Pkg. of 300
62797	RED, Pkg. of 12
62796	RED, Pkg. of 300
62670	ROYAL BLUE, Pkg. of 12
62680	ROYAL BLUE, Pkg. of 100
62690	ROYAL BLUE, Pkg. of 300
62655	PINK, Pkg. of 12
62664	PINK, Pkg. of 100
62675	PINK, Pkg. of 300

**.080" EVA Mouthguard Material Disc (120mm)**

Item #	Description
62814	CLEAR, Pkg. of 100
62812	CLEAR, Pkg. of 300
62671	ROYAL BLUE, Pkg. of 12
62681	ROYAL BLUE, Pkg. of 100
62691	ROYAL BLUE, Pkg. of 300
62666	PINK, Pkg. of 12
62676	PINK, Pkg. of 100
62686	PINK, Pkg. of 300
62791	RED, Pkg. of 300
62793	YELLOW, Pkg. of 300
62823	BLACK, Pkg. of 300

**.080" EVA Mouthguard Material Disc (125mm)**

Item #	Description
62807	CLEAR, Pkg. of 12
62808	CLEAR, Pkg. of 100
62809	CLEAR, Pkg. of 300
62813	WHITE, Pkg. of 300
62824	BLACK, Pkg. of 300
62835	ORANGE, Pkg. of 300

**.118" EVA Mouthguard Material**

Item #	Description
72530	WHITE, Pkg. of 12
72540	WHITE, Pkg. of 100
72550	WHITE, Pkg. of 300
72532	YELLOW, Pkg. of 12
72534	BLUE, Pkg. of 12
72536	ORANGE, Pkg. of 12
72539	ORANGE, Pkg. of 300
72538	RED, Pkg. of 12
72537	RED, Pkg. of 300
72500	BLACK, Pkg. of 12
72510	BLACK, Pkg. of 100
72520	BLACK, Pkg. of 300
62640	ROYAL BLUE, Pkg. of 12
62641	ROYAL BLUE, Pkg. of 100
62642	ROYAL BLUE, Pkg. of 300

**.118" EVA Mouthguard Material Discs (120mm)**

Item #	Description
62889	BLUE, Pkg. of 300
62888	CLEAR, Pkg. of 300
62890	RED, Pkg. of 300
62891	BLACK, Pkg. of 300
62899	WHITE, Pkg. of 300
62895	YELLOW, Pkg. of 100
62898	YELLOW, Pkg. of 300
62643	ROYAL BLUE, Pkg. of 300
62644	ROYAL BLUE, Pkg. of 300
62645	ROYAL BLUE, Pkg. of 300

**.118" EVA Mouthguard Material Discs (125mm)**

Item #	Description
62818	CLEAR, Pkg. of 12
62817	CLEAR, Pkg. of 100
62819	CLEAR, Pkg. of 300
62887	SKY-BLUE, Pkg. of 300
72630	BLACK, Pkg. of 12
72640	BLACK, Pkg. of 100
72650	BLACK, Pkg. of 300
72600	WHITE, Pkg. of 12
72610	WHITE, Pkg. of 100
72660	WHITE, Pkg. of 300
72620	BLUE, Pkg. of 300

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

[www.buffalodental.com](http://www.buffalodental.com)



**EVA Colored Mouthguard Sheets continued...**

**.160" EVA Mouthguard Material**

Item #	Description
62770	YELLOW, Pkg. of 12
62775	YELLOW, Pkg. of 25
62771	BLUE, Pkg. of 12
62777	BLUE, Pkg. of 25
62720	SKY-BLUE, Pkg. of 12
62740	SKY-BLUE, Pkg. of 100
62730	SKY-BLUE, Pkg. of 300
62772	RED, Pkg. of 12
62778	RED, Pkg. of 25
62773	ORANGE, Pkg. of 12
62779	ORANGE, Pkg. of 25
62774	BLACK, Pkg. of 12
62781	BLACK, Pkg. of 25
62783	BLACK, Pkg. of 300
62776	WHITE, Pkg. of 12
62782	WHITE, Pkg. of 25
62700	CLEAR, Pkg. of 12
62725	CLEAR, Pkg. of 25
62600	PINK, Pkg. of 12
62630	PINK, Pkg. of 100
62650	PINK, Pkg. of 300
626100	ROYAL BLUE, Pkg. of 12
626105	ROYAL BLUE, Pkg. of 100
626110	ROYAL BLUE, Pkg. of 300

**.160" EVA Mouthguard Material Disc (125mm)**

Item #	Description
62715	CLEAR, Pkg. of 12
62718	CLEAR, Pkg. of 100
62719	CLEAR, Pkg. of 300



MADE IN USA

**BUFF-TUFF Mouthguard Laminate**

Buff-Tuff Mouthguard Laminate features cross-graining for extra strength and constant thickness of finished mouthguards. Injection molded mouthguards cannot be cross grained! Contains *NO BPA, DIOXIN, PVC or INTERNAL BRACES!* Will not separate! **Designed Specifically For Athletic use.**

**BUFF-TUFF Mouthguard Laminate Squares & Discs**

Item #	Description
70000	BUFF-TUFF Clear 5"x5", Square, Pkg. of 10
70010	BUFF-TUFF Clear 5"x5", Square, Pkg. of 100
70100	BUFF-TUFF Clear 5"x5", Disc, Pkg. of 10
70110	BUFF-TUFF Clear 5"x5", Disc, Pkg. of 100
70120	BUFF-TUFF Blue 5"x5", Square, Pkg. of 10
70130	BUFF-TUFF Blue 5"x5", Square, Pkg. of 100
70140	BUFF-TUFF Blue 5"x5", Disc, Pkg. of 10
70150	BUFF-TUFF Blue 5"x5", Disc, Pkg. of 100
70160	BUFF-TUFF Orange 5"x5", Square, Pkg. of 10
70170	BUFF-TUFF Orange 5"x5", Square, Pkg. of 100
70180	BUFF-TUFF Yellow 5"x5", Square, Pkg. of 10
70190	BUFF-TUFF Yellow 5"x5", Square, Pkg. of 100
70200	BUFF-TUFF Black 5"x5", Square, Pkg. of 10



Item #	Description
70210	BUFF-TUFF Black 5"x5", Square, Pkg. of 100
70220	BUFF-TUFF Red 5"x5", Square, Pkg. of 10
70230	BUFF-TUFF Red 5"x5", Square, Pkg. of 100
70240	BUFF-TUFF Sky-Blue 5"x5", Square, Pkg. of 10
70250	BUFF-TUFF Sky-Blue 5"x5", Square, Pkg. of 100
70260	BUFF-TUFF White 5"x5", Square, Pkg. of 10
70270	BUFF-TUFF White 5"x5", Square, Pkg. of 100

**BUFF-TUFF Mouthguard Laminate Discs (125mm)**

Item #	Description	Item #	Description
70280	BUFF-TUFF Orange 5"x5", Disc, Pkg. of 10	70350	BUFF-TUFF Red 5"x5", Disc, Pkg. of 100
70290	BUFF-TUFF Orange 5"x5", Disc, Pkg. of 100	70360	BUFF-TUFF Sky-Blue 5"x5", Disc, Pkg. of 10
70300	BUFF-TUFF Yellow 5"x5", Disc, Pkg. of 10	70370	BUFF-TUFF Sky-Blue 5"x5", Disc, Pkg. of 100
70310	BUFF-TUFF Yellow 5"x5", Disc, Pkg. of 100	70380	BUFF-TUFF White 5"x5", Disc, Pkg. of 10
70320	BUFF-TUFF Black 5"x5", Disc, Pkg. of 10	70390	BUFF-TUFF White 5"x5", Disc, Pkg. of 100
70330	BUFF-TUFF Black 5"x5", Disc, Pkg. of 100	70400	BUFF-TUFF Clear 5"x5", Disc, Pkg. of 10
70340	BUFF-TUFF Red 5"x5", Disc, Pkg. of 10	70410	BUFF-TUFF Clear 5"x5", Disc, Pkg. of 100

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com



## Colored Mouthguard Materials with Polyshok™ Technology!

New Polyshok™ technology allows the creation of high-impact sports mouthguards that can absorb 150% more energy than standard EVA mouthguards. The improved shock absorption can result in less destructive forces being transferred to dental structures, fewer oral-fascial injuries, and less concussions.

### Item # Description

62901 .125" Polyshok™ Variety Pack, assorted box of 38 units

### Item # Description

62931 .155" Polyshok™ Variety Pack, assorted box of 30 units

### Item # Description

62902 .125" Black- Pkg. of 18

62903 .125" Blue- Pkg. of 18

62904 .125" Red- Pkg. of 18

62905 .125" Yellow- Pkg. of 18

62906 .125" Orange- Pkg. of 18

62907 .125" Purple- Pkg. of 18

62908 .125" Green- Pkg. of 18

62909 .125" White- Pkg. of 18

### Item # Description

62902 .125" Black- Pkg. of 18

62903 .125" Blue- Pkg. of 18

62904 .125" Red- Pkg. of 18

62905 .125" Yellow- Pkg. of 18

62906 .125" Orange- Pkg. of 18

62907 .125" Purple- Pkg. of 18

62908 .125" Green- Pkg. of 18

62909 .125" White- Pkg. of 18



## Polyshok™ GAMEguard (Pre-Formed "Boil & Bite)

High Performance Mouthguards w/ Carrying Cases- New! Polyshok™ GAMEguard SSD (No Strap) w/Case- Sports-Specific Design- Optimized Protection for Non-Helmeted Sports.

### Item # Description

62970 Large Red SSD Pre-Formed Mouthguard

62971 Large Blue SSD Pre-Formed Mouthguard

62972 Large Black SSD Pre-Formed Mouthguard

62973 Large Neon Green SSD Pre-Formed Mouthguard

62974 Large Neon Yellow SSD Pre-Formed Mouthguard

62975 Small Red SSD Pre-Formed Mouthguard

62976 Small Blue SSD Pre-Formed Mouthguard

62977 Small Black SSD Pre-Formed Mouthguard

62978 Small Neon Green SSD Pre-Formed Mouthguard

62979 Small Neon Yellow SSD Pre-Formed Mouthguard



Large Red SSD Pre-Formed Mouthguard

Polyshok™ GAMEguard PLUS w/ Strap & Case- Optimized Extra Protection for Helmeted Sports.

### Item # Description

62980 Large Red PLUS Pre-Formed Mouthguard

62981 Large Blue PLUS Pre-Formed Mouthguard

62982 Large Black PLUS Pre-Formed Mouthguard

62983 Large Green PLUS Pre-Formed Mouthguard

62984 Large Yellow PLUS Pre-Formed Mouthguard

62985 Small Red PLUS Pre-Formed Mouthguard

62986 Small Blue PLUS Pre-Formed Mouthguard

62987 Small Black PLUS Pre-Formed Mouthguard

62988 Small Green PLUS Pre-Formed Mouthguard

62989 Small Yellow PLUS Pre-Formed Mouthguard



Large Red PLUS Pre-Formed Mouthguard

### Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA

Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610

www.buffalodental.com