

# Opti™ Workstation

The ideal lab workstation with built-in dust collection



## Powerful suction in a small size!

The Opti™ Workstation is a brightly illuminated, easy to use, small dust collecting cabinet for use with the Optiblast II, electric handpieces, or other dust-producing devices in the office or the lab.

Its small size makes it suitable for bench tops and office back rooms.

### Features and Benefits:

- Small size (11" W x 12" D x 8 1/2" H) fits easily on table tops
- Hinged top window and internal light fixture provide bright illumination for clear viewing
- Dual, wide hand holes allow comfortable hand positioning and access to large workspaces
- Efficient suction fan and filter collect damaging grinding debris and maintain quiet operation
- Additional 120 V AC power outlet on the back of the unit for operating electric devices
- Hinged rear door for easy cleaning



### Ordering Information:

<b>Item#</b>	<b>Description</b>
36560	Opti Workstation w/ Suction, 120 V AC
36570	Opti Workstation w/o Suction, 120 V AC (w/ 2 7/8" OD port on back for connection)
36565	Opti Workstation w/ Suction, 220 V AC
36575	Opti Workstation w/o Suction, 220 V AC (w/ 2 7/8" OD port on back for connection)

### Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203

Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com